(DMH6000)

OZONE HANDY MISTER

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The background features a light blue gradient with several 3D DNA double helix structures on the left side, rendered in white and light blue. On the right side, there are glowing blue network graphics consisting of interconnected nodes and lines, resembling a molecular or data network. The overall aesthetic is clean, modern, and scientific.

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Sterilization Range and Effect: 滅菌範圍と効果

OZONR HANDY MISTER

Certificate and Test Report

Scan the QR code
Check the original
inspection report

OZONR HANDY MISTER

Certificate and Test Report

Microbiological Test Results		Elimination Rate
Microbiological Testing	Cutibacterium acnes	99.99%
	Staphylococcus epidermidis	99.99%
	Staphylococcus aureus	99.94%
	Streptococcus	99.99%
	Malassezia	99.99%
	Pseudomonas aeruginosa	99.95%
	Escherichia coli	99.92%
	Mite Removal Rate	99.00%

*Data sourced from tests conducted by SGS-CSTC Standards

Technical Services (Ningbo) Co., Ltd. and Pony Testing.

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Four Core Advantages

Multi-Functionality & Integrated Skincare 多機能性&統合スキンケア

Beyond basic hydration, this device integrates multiple skincare functions including sterilization, anti-inflammation, skin repair, and antioxidant protection.

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Four major cores

Safe & Zero Additives Naturally Gentle 安全で無添加自然な優しさ

It uses only water, with no alcohol, fragrance, preservatives or other chemical additives. Non-irritating and side-effect-free, it is suitable for all skin types, including sensitive skin.

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
On-Demand Generation Eco-Friendly & Cost-Effective

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3.2

Comparison with Market Competitors

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Business Positioning and Value

Positioning: 位置付け

It is far more than a home-use beauty device; it is a water-centric at-home health and skincare solution, representing a pivotal extension of EOG technology from "cleaning" to "skincare care" scenarios.

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Value: 価値

Fills a critical market gap — positioned between entry-level mist sprayers (priced at hundreds of RMB) and professional high-end oxygen infusion devices (priced at thousands of RMB), it delivers a truly bioactive, safe, and cost-effective intelligent skincare choice.

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Women aged 25-40 who pursue safe, highly effective skincare and are tired of excessive chemical ingredients in their routines

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Individuals seeking a physical approach to inhibit bacteria at the source and regulate the skin's water-oil balance.

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Those in need of a one-stop, fast and effective daily skincare routine.

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5

Golden Sales Scripts & Practical Strategies

5.1

Pain Point-Focused Opening Script

Scenario: For customers with acne-prone skin, sensitive skin, or those concerned about skincare product ingredients.

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"Have you ever experienced this: your skin becomes sensitive and red with seasonal changes, many skincare products don't absorb well, and you're prone to clogged breakouts? Or do you worry that complex skincare ingredients are burdening your skin?"

Today, I'm not recommending a new bottle of serum — I'm introducing an Active Oxygen Hydration Device that lets you "skincare with tap water". It uses zero chemical additives. Relying on the natural active oxygen generated by electrolyzing water, it helps you sterilize, reduce inflammation at the source, and deeply hydrate your skin, bringing it back to a healthy, balanced state."

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User Pain Points	Competitor Limitations	Our Solution (Sales Script)
<p>"My skin feels dry right after using a hydrating mist "</p> <p>「保湿ミストを使った直後、肌が乾燥する感じがする」</p>	<p>Traditional mist sprayers produce large water droplets that only moisten the stratum corneum, leaving skin even drier after evaporation.</p> <p>従来のミストスプレーは大きな水滴を発生させるため、角質層を湿らせるだけで、蒸発後には肌がさらに乾燥してしまう。</p>	<p>"Our nano-mist is just 0.3 microns, 1/200 the thickness of a human hair. It sprays directly into the dermis to fully hydrate it, not just 'coat the surface' with water. Remember to lock in moisture after use, and the effect lasts all day long. "</p> <p>当社のナノミストはわずか0.3ミクロン、人間の髪の毛の200分の1の細さ。表面を水で覆うのではなく、真皮に直接噴射してしっかり保湿。使用后、保湿剤でしっかり保湿すれば、効果は一日中持続します。</p>
<p>"My acne keeps coming back, and I'm afraid to use random products "</p> <p>ニキビが何度も再発するので、適当な製品を使うのが怖いんです。</p>	<p>Many anti-acne products contain acids and alcohol, are highly irritating, and only treat symptoms, not the root cause.</p> <p>多くのニキビ治療薬には酸やアルコールが含まれており、刺激が強く、症状を緩和するだけで根本原因には対処できません。</p>	<p>"Acne is fundamentally caused by bacteria. The active oxygen this device emits is your skin's 'natural anti-inflammatory water'. It gives your pores a sterile clean every day, inhibiting bacteria at the source so acne has no chance to come back. "</p> <p>「ニキビの根本原因は細菌。この機器が放出する活性酸素は、『天然の抗炎症水』です。毎日毛穴を清潔に保ち、細菌の発生を抑制し、ニキビの再発を防ぎます。」</p>
<p>"I have sensitive skin and am allergic to almost everything "</p> <p>「私は敏感肌で、ほとんどすべてのものにアレルギーがあります。」</p>	<p>Most products on the market have complex ingredients, requiring careful selection.</p> <p>市販されている製品のほとんどは複雑な成分を含んでいるため、慎重な選択が必要となる。</p>	<p>"Its biggest advantage is simplicity: it only uses water. After electrolysis, the active oxygen breaks down into oxygen and water in minutes, leaving zero chemical residue. Many beauty salons in Japan even use ozone treatments for sensitive clients, so the safety level is extremely high. "</p> <p>「最大の利点はそのシンプルさ。使用は水だけ。電気分解後、活性酸素は数分で酸素と水に分解され、化学物質の残留物は一切残りません。日本では多くの美容サロンが敏感肌のお客様向けにオゾン療法を取り入れ、安全性は非常に高いと言えます。」</p>
<p>"Beauty devices are too expensive, and I have to buy special serums too "</p> <p>美容機器は高すぎるし、専用の美容液も買わなきゃいけないし。</p>	<p>High-end oxygen infusion devices are pricy, and the ongoing cost of purchasing serums is enormous.</p> <p>高性能の酸素注入装置は高価であり、血清を購入し続ける費用も莫大だ。</p>	<p>"Our device costs just ¥399, and after that, you only need tap water. It doesn't rely on expensive serums; it generates its own active ingredients. Do the math: you'll pay for the device with the money you save on just one bottle of serum. This is smart skincare. "</p> <p>「このデバイスは¥で、あとは水道水だけ。高価な美容液に頼る必要はなく、有効成分を自ら生成します。美容液1本分の節約でこのデバイスの費用を回収できます。」</p>

The background features a light blue gradient with several 3D DNA double helix structures on the left side and abstract network graphics consisting of glowing blue nodes and lines on the right side.

6

Usage and Maintenance Instructions

6.1

Operation Steps

HOW TO USE

1

STEP1 Open the water tank

2

STEP2 Add water

3

STEP3 Twist the water tank

4

STEP4 Turn on the sprayer

6.2 Daily Maintenance

Regular cleaning and maintenance with citric acid every 3-4 months: Mix citric acid with warm water at a ratio of 1:5 (e.g., 2 grams of citric acid with 10 milliliters of warm water at approximately 40°C). After the citric acid is fully dissolved, soak the module with the solution for 5-10 minutes. Then pour out the citric acid solution and rinse the module thoroughly with clean water.

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- Regular cleaning of the water tank: Twist the water tank out clockwise and clean the inner wall of the water tank with a soft brush or sponge to remove limescale and dirt.

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- Long-term storage: If the device is not in use for an extended period, clean up limescale and dirt first and store it in a dry environment.

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